



OSSEO UNO  
SHAVING SYSTEM

# OSSEOUNO, TO MEET THE MOST DEMANDING REQUIREMENTS.

---

The new Bien-Air shaver control offers you cutting-edge features and small size. Fitted with the S120 handpiece, the system is easy to use, has a high performance level and the legendary reliability that is a feature of all Bien-Air products.

The **OSSEOUNO** system has been designed for shaping bones and for resection of soft and hard tissues as part of eye, nose and throat surgical operations.

All the functions of the **OSSEOUNO** control system can be viewed and accessed by simply pressing keys or the pedal.

The **OSSEOUNO** motor can be programmed in continuous rotation or oscillating rotation mode. You have a choice of five speeds.

The **OSSEOUNO** is provided with an integral peristaltic pump. Pump delivery is easily adjustable from 15 to 130 ml/min. by simply pressing a key.

## COMFORT AND PRECISION

The S120 handpiece, with integral micro-motor, is characterised by its lightness, precision and handling comfort. Its unique ergonomics and its integral irrigation system ensure perfect ease of handling. Even in hard-to-access areas, it offers high precision operation. The S120 is provided with a straight suction channel. This channel allows continuous operation without clogging.

The S120 is coupled to a brushless motor without sensors which is renowned for its power and reliability. The S120 allows blade to be positioned at any desired angle while operating. Meeting the most stringent hygiene standards, the handpiece and its cable can be sterilised in an autoclave in complete security.

The S120 handpiece allows various types of disposable or reusable shaver blades to be connected using just one hand.

OFFER YOURSELF  
SWISS QUALITY.

With 50 years' experience in ENT surgery and dentistry, Bien-Air offers a range of solutions and products which are renowned worldwide for their excellent quality.

100% Swiss Made, Bien-Air products are designed and manufactured in the heart of Watch Valley, a region famed for its unique expertise in microtechnology. Our products have been developed directly from the daily work of professionals practising in the surgical field.

Our customers testify daily to the high quality, reliability and precision of our products, on terms that are more attractive than ever.



The OSSEOUNO system will convince you through its great ease of use, efficiency and reliability.

Bien-Air Surgery SA  
Rue de l'Ouest 2b  
CH-2340 Le Noirmont  
Switzerland  
Tél. +41(0)32 344 64 40  
Fax +41(0)32 344 64 45  
surgery@bienair.com

Bien-Air  
Deutschland GmbH  
Surgery  
Jechtinger Strasse 11  
79111 Freiburg, Deutschland  
Tel. +49 (0)761 45 57 40  
Fax +49 (0)761 47 47 28  
ba-d@bienair.com

Bien-Air España SA  
Surgery  
Entença, 169 Bajos  
08029 Barcelona, España  
Tel. (+34) 934 25 30 40  
Fax (+34) 934 23 98 60  
ba-e@bienair.com

Bien-Air USA, Inc.  
Medical Technologies  
5 Corporate Park  
Suite 160  
Irvine, CA 92606, USA  
Phone 1-800-433-BIEN  
Phone 949-477-6050  
Fax 949-477-6051  
ba-usa@bienair.com

Bien-Air France Sarl  
Surgery  
55-57, avenue Jean Lolive  
93500 Pantin Cedex, France  
Tel. +33 (0)1 41 83 60 70  
Fax +33 (0)1 48 96 07 40  
ba-f@bienair.com

Bien-Air on Internet:  
www.bienairsurgery.com

Bien-Air Italia s.r.l.  
Surgery  
Via Vaina 3  
20122 Milano, Italia  
Tel. +39 (02) 58 32 12 51/52/54  
Fax +39 (02) 58 32 12 53  
ba-i@bienair.com

Bien-Air UK Limited  
Arundel House / Unit 1 - Ground floor  
Amberley Court / Whitworth Road  
Crawley, West Sussex,  
RH11 7XL, England  
Tel. +44 (0)1293 550200  
Fax +44 (0)1293 520481  
ba-uk@bienair.com

Bien-Air Asia Ltd.  
Surgery  
Nishi-Ikebukuro  
Daiichi-Seimei Bldg. 10F  
2-40-12 Ikebukuro, Toshimaku  
Tokyo, 171-0014, Japan  
ビエン・エア・アジア株式会社  
〒171-0014  
東京都豊島区池袋2-40-12  
西池袋第一生命ビルディング10F  
Tel. +81 (3) 5954-7661  
Fax +81 (3) 5954-7660  
ba-asia@bienair.com

Beijing Bien-Air  
Medical Instrument  
Technology Service Co. Ltd.  
Room 907, The Exchange Beijing,  
No 118 Jian Guo Lu Yi,  
Chao Yang District,  
Beijing 100022, China  
北京彼岸医疗器械  
技术服务有限公司  
北京市朝阳区建国路  
乙118号招商局中心  
京汇大厦907室  
电话: +86 10 6567 0651  
传真: +86 10 6567 8047  
surgery-beijing@bienair.com

## Technical data



### OSSEOUNO

Marking	CE 0120 as per 93/42 CEE directive
Supply voltage	100 - 240 Vac / 2A - 1A / 50/60 Hz
Motor power supply	Quick connector of the push-pull type / 48 Vdc
Electrical insulation class	Class I
Parts applied	Type BF
Medical class	Ila
Protection class	IP 41
Dimensions / Weight	(L x H x D) 261 x 110 x 300 mm (height with irrigation stand 520 mm) / 4,1 kg
Peristaltic pump	Delivery of 15 to 130 ml/min., adjustable by 30% steps (with two-speed Bien-Air Surgery irrigation line REF 1100037-001)
Irrigation stand	Stainless steel, height 450 mm



### S120 handpiece with cable (REF 1700302-001)

Type	Direct current, without brushes or Hall detectors
Torque	0 to $\pm 15$ mNm
Speed in continuous rotation	500 to 8000 rpm
Speed in oscillation mode	500 to 5000 rpm
Size / Weight	Diameter 18 mm, length 140 mm, angle 15 / 198 g without cable
Motor cable	Push-pull type at <b>OSSEOUNO</b> end, length 295 cm $\pm 5$ cm
Instruments	Straight and angled shaver blades, disposable or reusable
Sterilisation	In autoclave at 134°C for 6 min.



### Multifunction pedal

Protection class	IPX8
Functions	Selection of 5 speeds in motor clockwise or oscillating rotation, variable speed drive, peristaltic pump start/stop and selection of 3 delivery rates
Dimensions / Weight	(L x H x D) 160 x 55 x 170 mm / 0,830 kg
Cable	295 cm $\pm 5$ cm

**OSSEOUNO set**

**Shaver S120**

**REF 1700300-001**

The company reserves the right to make model modifications and technical changes.